BioResearch

LONZA

Nucleofector™ Technology

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Nucleofector™ Technology
for Hard-to-transfect Cells

Introduction: Nucleofector™ Technology

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Nucleofector™ Technology – the Superior Non-viral Method


Figure 1. Transfection of the human natural killer cell line NKL using traditional electroporation and Nucleofection. 5×10^5 NKL cells were transfected with 2.5 μ g of pmaxGFP™ Vector. Nucleofection: Nucleofector™ Solution V; Program O-017. Competitor B electroporation: 25 mV, 96 μ F. Transfection efficiency was monitored by flow cytometry after 24 hours. Cells transfected by Nucleofection show a significantly better transfection efficiency compared to cells transfected by traditional electroporation. Cell viability, as measured 18 hours after transfection, was also superior using Nucleofection. (Data courtesy of Dr. John Coligan, Laboratory of Immunogenetics, NIH/NIAD, Rockville, MD, USA. J Immunol Methods [2004] 284: 133-140.)

DNA Delivery Straight Into the Nucleus



Figure 2. Normal human dermal fibroblasts (neonatal) were transfected with 2.5 μ g TMR-labeled plasmid DNA encoding eGFP. After 2 hours, cells were fixed with 3.5% PFA and analyzed by confocal microscopy. TMR label is shown in (A), GFP fluorescence in (B), DAPI nuclear staining in (C) and a merge of all three fluorescent labels in (D).

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Overview About Nucleofection Platforms

	Advanced Platform	96-well Add-on	High-throughput Platform	Basic Device
Device	4D-Nucleofector™ System	96-well Shuttle® Device	HT Nucleofector™ System	Nucleofector™ 2b Device
				
Unit				
Throughput (samples per run)	Low to medium (1-16)	Low to high (1-96)	High (384)	Low (1)
Reaction volume	20 µL + 100 µL	20 µL	20 µL	100 µL
Electrode material	Conductive polymer	Conductive polymer	Conductive polymer	Aluminum
Low cell numbers (20 µL)	2×10^4 to 1×10^6	2×10^4 to 1×10^6	2×10^4 to 1×10^6	—
High cell numbers (100 µL)	2×10^5 to 2×10^7	—	—	2×10^5 to 2×10^7
DNA Vector amount/sample	0.2 – 1 µg (20 µL) 1 – 5 µg (100 µL)	0.2 – 1 µg	0.2 – 1 µg	1 – 5 µg
siRNA amount/sample (concentration 2nM–2µM)	0.2 – 200 pmol (100 µL) 0.04 – 40 pmol (20 µL)	0.04 – 40 pmol	0.04 – 40 pmol	0.2 – 200 pmol
Adherent Nucleofection (20 µL Nucleocuvette™ Strips)			—	—
Adherent Nucleofection (24-well culture plates)		—	—	—
Compatibility with 96-well Shuttle® Device		—	—	—

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Transferability Between Nucleofection Conditions Between Different Formats


Figure 6. Various primary cells were transfected in the two Nucleocuvette™ vessel formats (20 µL and 100 µL) using the indicated programs. Twenty-four hours post Nucleofection cells were analyzed for transfection efficiency (flow cytometry) and viability (cell number normalized to no program control).

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
Figure 7. Efficient adherent Nucleofection of neurons in 24-well culture plates. Mouse cortical neurons were seeded into poly-D-lysine coated 24-well plates (1×10^5 cells/well). After 6 DIV, cells were transfected with pmaxGFP™ Vector using the AD1 4D-Nucleofector™ Y Kit. One day post Nucleofection, cells were stained by MAP2 antibody (red) and analyzed by fluorescence microscopy for maxGFP™ protein expression.



Figure 8. Human umbilical vein endothelial cells (HUVEC) were isolated and plated in passage 1 into collagen-coated 24-well plates at a density of 50,000 cells/well. After 1 DIV cells were transfected with 16 µg pmaxGFP™ Vector using AD1 4D-Nucleofector™ Y Solution and program CA-215. Cells were analyzed for maxGFP™ Protein expression after 24h. (Data kindly provided by M. Sauvage, Pharmaceutical Industry, FR)

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High Transfection Efficiencies in Suspension Cell Lines


Figure 11. Transfection efficiencies 24 hours post Nucleofection in selected cell lines relevant for immunology research. Cells were transfected with either eGFP, maxGFP[®] Reporter Protein or H-2KK and analyzed 24 hours post Nucleofection. Viability ranged from 60–90%.

Ordering Information

Description	Cat. No.	
Nucleofector™ Devices		
4D-Nucleofector™ Core Unit	AAF-1001B	—
4D-Nucleofector™ X Unit	AAF-1001X	Requires core unit to build complete system
4D-Nucleofector™ Y Unit	AAF-1001Y	Requires core unit to build complete system
4D-Nucleofector™ Guarantee, 2-year extension	AWE-1002	Has to be purchased at the time the device is purchased
4D-Nucleofector™ Service Contract	AWC-1001	Can be purchased at any time during the guarantee period
96-well Shutte™ add-on (including laptop)	AAM-1001S	Requires core and X unit to build complete system
96-well Shutte™ Guarantee, 2-year extension	AWM-1002	Has to be purchased at the time the device is purchased
96-well Shutte™ Service Contract	AWB-1001	Can be purchased at any time during the guarantee period
HT Nucleofector™ System	AAU-1001	—
HT Nucleofector™ Installation and Training	AWT-1001	—
HT Nucleofector™ Service Contract	AWU-1001	Can be purchased at any time during the guarantee period
Nucleofector™ 2b Device	AAB-1001	—
Nucleofector™ 2b Guarantee, 2-year extension	AWD-2002	Has to be purchased at the time the device is purchased
Nucleofector™ 2b Service Contract	AWA-2001	Can be purchased at any time during the guarantee period

	100 µL Nucleocuvette™		20 µL Nucleocuvette™; 16-well	Dipping Electrode
	12 rxn	24 rxn	32 rxn	24 rxn
4D-Nucleofector™ Kits				
P1 Primary Cell 4D-Nucleofector™ X Kit	V4XP-1012	V4XP-1024	V4XP-1032	—
P2 Primary Cell 4D-Nucleofector™ X Kit	V4XP-2012	V4XP-2024	V4XP-2032	—
P3 Primary Cell 4D-Nucleofector™ X Kit	V4XP-3012	V4XP-3024	V4XP-3032	—
P4 Primary Cell 4D-Nucleofector™ X Kit	V4XP-4012	V4XP-4024	V4XP-4032	—
P5 Primary Cell 4D-Nucleofector™ X Kit	V4XP-5012	V4XP-5024	V4XP-5032	—
Primary Cell Optimization 4D-Nucleofector™ X Kit	—	—	V4XP-9096 (96 rxn)	—
Basic Neuron 4D-Nucleofector™ X AD Kit (adherent)	—	—	V4XP-1A32	—
AD1 4D-Nucleofector™ Y Kit (adherent)	—	—	—	V4YP-1A24
AD2 4D-Nucleofector™ Y Kit (adherent)	—	—	—	V4YP-2A24
Optimization 4D-Nucleofector™ Y Kit (adherent)	—	—	—	V4YP-9A48
SE Cell line 4D-Nucleofector™ X Kit	V4XC-1012	V4XC-1024	V4XC-1032	—
SF Cell line 4D-Nucleofector™ X Kit	V4XC-2012	V4XC-2024	V4XC-2032	—
SG Cell line 4D-Nucleofector™ X Kit	V4XC-3012	V4XC-3024	V4XC-3032	—
Cell Line Optimization 4D-Nucleofector™ X Kit	—	—	V4XC-9096 (64 rxn)	—

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